



Infant Protection Solutions

Security that evolves with your organization's needs

KinderGUARD® infant protection solutions serve as a personal guard for each infant in your care. The small ankle transmitter is a visible sign to parents and loved ones of your commitment to safety and security. KinderGUARD® offers peace of mind for new parents and staff while elevating your facility's performance and reputation.

KinderGUARD® is specifically designed to grow alongside your facility's evolving needs, empowering you with a higher degree of flexibility and adaptability. With its user-friendly interface, KinderGUARD® is easy to manage and helps protect your most precious patients at all times.



Universal Series Features

Dependable by Design

KinderGUARD Universal series is computer-enhanced, not computer-dependent – so your tiniest patients are protected, even in the event of power outage or computer downtime

Automatic Enrollment

Once applied to a newborn or infant, a transmitter automatically becomes active in the system ensuring immediate protection

Total Supervision

Active transmitters report into the system every 7 to 10 seconds to insure supervision around the clock

Newborn and Infant Friendly

Our straps are made of an advanced ultrasoft material for the highest level of comfort against delicate skin

Tamperproof and Accurate

Leading-edge cutband technology sends an alert to staff if a strap is cut, broken or compromised. A design superior to traditional skin sensing tags which can send false alerts due to loose fit or quick body shift



The Right Technology to Protect the Most Vulnerable

A KinderGUARD® infant protection system uses radio frequency (RF) technology which offers greater reliability and security than one relying on Wi-Fi. KinderGUARD® operates on dedicated frequencies that are less prone to network interference, ensuring continuous monitoring even if the hospital's Wi-Fi network experiences outages or congestion.

KinderGUARD® solutions don't depend on internet connectivity, reducing latency and cybersecurity risks. The systems provide more precise real-time location tracking within controlled zones, KinderGUARD® ideal for safeguarding infants in critical hospital environments.



1125



Universal Series Features *(cont.)*

Ease of Use

System provides quick reference visuals of alarm/alert details, expired transmitter notifications and a clear display of monitored floorplan

Integration Friendly

Seamless integration with card access control, nurse call, CCTV and more

Secure Performance

All KinderGUARD® tags have a unique ID and operate on select frequencies to protect against RF interference

Efficient and Ergonomic

KinderGUARD® tags are unobtrusive, light weight, easily sanitized and reassignable throughout their lifespan.



System Software Enhanced Features

- Quick, straightforward, and easy to use
- Near real-time status
- Virtual environment supported
- Clear visual display of floor plans and exits
- Simple infant transfer
- Multiple-floor monitoring capabilities
- Displays detailed alarm and alert information including:
 - Transmitter ID
 - Event location
 - Infant name
 - Room number
- Quick report generation including alarm/alert histories, resident activity logs, staff response times, expired transmitter list, transmitter re-order information



Secure Care Products

Helping you enhance the safety of your patients and staff with our line of selective monitoring solutions