

# Universal Mini-Server for Wander Management

Integrated software solution unlocking next-level safety capabilities

The **cost-conscious Universal Mini-Server hardware solution** allows facilities with small to mid-size wander management needs to **enable the full range of enhancements offered by the DoorGUARDIAN® ID system software.**

**Universal Mini-Servers are self-contained** and can eliminate the need for client sourced server and LAN network solutions, and are **scalable to accommodate each client's unique needs.**



## ID System Software Features

- Intuitive, easy to use onscreen interface
- Near real-time status
- Clear visual display of floor plans and exits
- Simple resident transfer
- Multiple-floor monitoring capabilities
- Displays detailed alarm and alert information including:
  - Tag ID
  - Event location
  - Resident name
  - Room number



- Quick report generation including:
  - Alarm/alert histories
  - Expired tag list
  - Tag activity logs
  - Tag reorder information
  - Staff response times



## Universal Mini-Server Components

**Self-contained, turnkey and scalable**



Mini  
Server



Wi-Fi  
Router



EXIU  
Hub

*A reliable, plug-and-play solution that enables the capabilities of DoorGUARDIAN® ID System Software with **minimal budget impact.***

*Perfectly suited for small to mid-size facilities, **Universal Mini-Servers can be scaled to address the needs of multiple locations within a facility** and accommodate a variety custom communication/reporting requirements.*

1125



## ID System Software Practical Applications *(a few of many)*

**Tracking the movement of wander-prone residents** helps staff quickly detect when someone is at risk of leaving safe areas, reducing the likelihood of injury or getting lost.

The location data recorded also allows **caregivers to respond faster in emergencies**, improving overall resident safety and well-being along with greater peace of mind for family members.

Additionally, movement data can highlight resident behavior patterns, **enabling personalized care plans** that better support each resident's needs.

Monitoring the location of staff within the facility helps **ensure that caregivers are efficiently distributed**, reducing response times and improving resident care.

The information provided also **enhances accountability and transparency**, making it easier to verify that critical tasks and rounds are completed.

Additionally, analyzing staff movement patterns can **identify workflow bottlenecks**, leading to better scheduling and more effective use of resources.



## Universal Mini-Server Specs at a Glance



Mini Server

- 3.4 GHz Celeron processor
- 32 GB DDR4 RAM
- 1 TB SSD Hard Drive
- Intel Graphics 1.2 GHz
- 32 GB Graphics Card RAM
- 5.5 x 5.5 x 1.5 inches
- 1.1 lbs
- Supports multiple device connections



Wi-Fi Router

- Ultra-fast Wi-Fi® 6 technology
- 802.11ax
- Wi-Fi Security
- Four Gigabit ethernet ports
- One Gigabit internet port
- One USB 3.0 port- ext. HD support
- Wall mountable
- 8.5 x 6.1 x 2.2 inches
- 1.5 lbs



EXIU Communication Hub

- Manages and passes information from DoorGAURDIAN® hardware to the PC-based software to initiate appropriate actions and alerts.
- 9.0 x 5.5 x 1.0 inches
  - Surface or flush mount
  - 12V DC or 1 Amp Plug-in adapter
  - Onboard self-tests
  - Battery Back-up