



# Infant Protection Solutions



## Automatic Enrollment

When a transmitter is applied to a child, it automatically becomes active in the system offering instant protection.

## Total Supervision

Transmitters are active and report in every 7 to 10 seconds.

## Reusable Transmitters

Our transmitters are small, light-weight, hypoallergenic and easy to clean for reuse.

## Tamper-Proof

Straps feature cutband technology that will alert staff if cut, broken or tampered with.

## Baby and Child Friendly

The adjustable strap is soft and comfortable against delicate skin.

## Easy-To-Use Software

Patient management simplified, see alarm/alert details, expired transmitter notifications and clear display of monitored floorplan

## Reliable Security

KinderGUARD is not computer controlled or dependent, so it continues providing protection even in the event of a computer system failure or power outage.

## Technical Support 24/7/365

Get support anytime from your local certified distributor and Secure Care Products' customer support line.



# Bring Your Facility Into the Future of Infant and Pediatric Care With



Real-Time Locating



Mobile Alerting

## Be Connected Anywhere You Go


Adding the mVision mobile app to your solution enables customized communications to your staff with real-time notifications that include picture push directly on their iOS or Android mobile devices.

## Quickly Locate Key Personnel & Assets

Be able to find the location of all monitored people and assets within seconds. Reduce time wasted manually searching the unit for whatever or whoever you need.

## Seamlessly Integrate with KinderGUARD Infant Protection

Real-time location information enables you to use data to review and improve how you protect your patients, staff, assets and automate workflows.



# Secure Care

Products, LLC

For more information on our solutions visit:

[www.securecare.com](http://www.securecare.com)

39 Chenell Drive, Concord, NH 03301

[info@securecare.com](mailto:info@securecare.com)

1 (800) 451-7917 or (603) 223-0745



Made in the USA

Meets UL Standards

