


# ID Based Wander Management Solutions



## Improve Staff Response Times

A Secure Care ID-based system allows you to link resident or patient information to the transmitter they wear. This means if an alarm or alert sounds, you can quickly identify who they are and determine the exact number of people involved.

The Secure Care software will display a floor map letting you know exactly which area of the facility an event is happening in. You can also associate a picture with each transmitter during registration. That way when a person is involved in an event, the software will display a visual of the resident/patient involved, enhancing recognition and response time.

### Get Information sent directly to mobile devices

mVision® allows staff members to receive notifications of events as they occur. This way they do not need to be near the software or an exit to hear an alarm occur. They will be notified with all the pertinent information sent directly to their smartphone or tablet, leading to faster and more targeted response to emergencies.



Wandering  
Residents

Wandering  
Patient

Location-Based  
Data Solutions

Senior Care  
LBDS

Enterprise  
Solutions

### About Secure Care Products

Global leader in the design, manufacture, distribution, and service of specialty RFID products providing security solutions to healthcare and senior care facilities

Provide solutions for infant security, wandering resident/patient, mobile communication and location-based data solutions.

Products featured in hospitals, long-term-care facilities, assisted living communities, memory-care and sub-acute facilities.

Exceptional service and support to our customers and partners.

Highest quality products in order to achieve the safest environment for patients/residents.

Ease of use for our end-users, and instill peace of mind for facility staff and patient/resident family members.

For more information on our solutions visit:

[www.securecare.com](http://www.securecare.com)



A Subsidiary of Valstone Systems Corporation USA, Inc.